

PCEDX Navi

File View Tools Maintenance Help

Navigation Panel

- Start-up
- Initialization
- Cooling Detector
- X-ray ON
- Stabilizing
- Check Analysis
- Energy Check
- Control Analysis
- Analysis
- Preparation
- Measurement
- Results
- Shutdown

## Results : Connector

PCEDX Navi

Element	Judgment	Result	Unit	
Cd	OK	ND	ppm	14.2
Pb	OK	ND	ppm	7.7
Cr	OK	52.2	ppm	10.7
Hg	OK	ND	ppm	9.0
Br	NG	36727.0	ppm	452.7
Cl	OK	ND	ppm	60.9
Sb	NG	10696.9	ppm	211.3
Plastic		95.252	%	0.000

Sample Information

Sample Name : Connector


Meas. Date : 2021-06-29 13:00:43

Group : [ Qual-Quant. ] ScreeningFast

Comment: Fast mode with Material check


Operator : Shimadzu

Sample Image :




Next Sample Result List Report..

Ready Ready 50KV 30uA Filter: DT%30 Collimator:10mm



Comparison of Lead Profiles in Copper Alloys



Comparison of Chromium Profiles in Copper Alloys